

EXHIBIT "A"

11/13/01 -11/27/01

EB on Scaffold Expt: Comparison of RA, TGFB, and Actavin Growth Factors

METHODS:

Preparation of Scaffolds:

1. Sterilized 10 polymer scaffolds overnight in 70% ETOH.
2. Soaked scaffolds in 3 changes of PBS, about 5 mins each change. Transferred with sterilized forceps.

Preparation of EBs:

3. Transfer EBs in media to 50 mL Falcon tube. Pipetted lightly to remove EBs adhered to plate. Washed plate with 5 mL old media.
4. Removed most supernatant. Resuspended in 5 mL new EB media. Transferred to 15 mL Falcon tube.
5. Centrifuged at 800 for 1 minute with brake.
6. Diluted mef trypsin 1:5 (white: 25g porcine tryp/L in 0.9% NaCl) in PBS. Add 2 mL diluted trypsin to the 10mL EB/media mixture. Resuspended by pipetting.
7. Incubated for 5 mins. Added 4mL TNS to cells. Resuspended, and centrifuged at 800 for 3 min. Remove supernatant. Add EB media to aliquot the cells into 2:1 (for scaffolds/EB: EB samples)

Scaffold Conditions:

GF	Stock	Working	# scaffolds	EB Plates w/o Scaffold
RA	7.5mg/mL in DMSO	(1:25000)	2	1
TGF-B	2ug/mL	(1:1000)	2	1
Activin	20ug/mL	(1:1000)	2	1
Activin & TGF-B		(1:1000)	2	1
Normal media			2	1

8. Prepared 25mL EB media with GFs. (RA 1:25000; TGF-B 1:1000; Act 1:1000)
9. Added 1mL of respective media to each well to be used for scaffolds. Added scaffolds to wells.
10. Prepared matrigel 1:1 with respective media (50 uL matrigel with each media .
11. Aliquotted cells into eppendorfs. Spin down <1000 for 4 min. Remove supernatant carefully. Add 20-25 uL matrigel-media to cells. Mix well.
12. Removed media from scaffolds in wells.
13. Added 20-25 uL cell-matrigel mixture onto each scaffold.
14. Incubate 30 mins to solidify matrigel.
15. Added 4 mL respective media to each well.
16. For remaining EB's aliquot, spin down and add media. Pipette into wells and add respective media.
17. Placed on shaker.

10

11/3

EB on Scaffold with GF Expt

Scaffold - PLG-A

70%

1. Sterilize polyimide scaffold in EtOH overnight
2. Rinse in PBS: Pour out until only 10ml EtOH. Pour PBS in small wells. Pour all scaffold + EtOH into plate. Transfer into PBS in 6-well plates with flickup.
3. Transfer 3k per ~5min each

EB (pyelectro)

4. Transfer EB + media to 50ml Falcon. Pipette lightly to remove EBs stuck to plate. Wash with clean media or old media.
5. Remove supernatant until 5ml left. Resuspend in new media (10ml). Transfer to 15ml Falcon. Use same 5ml for all 5 tubes.
6. Repeat again with 5ml. S. Centrifuge 800 for 1min with brake.
7. Dilute Trypsin (white, 10X): 1 → 5. Add 2ml diluted to the 10ml EB-media mixture. Resuspend by pipetting. Pour into 2 small petri dishes (5ml each).
8. Incubate 3min in the 15ml Falcon

[Dilute Trypsin 1:5 using PBS → only for EB's]

Scaffold Conditions = 10

1. RA ~~10~~ 2
2. TGFβ ~~10~~ 2
3. Acturin ~~10~~ 2
4. TGFβ + Acturin ~~10~~ 2
5. control 2

EB's w/o scaffold

- TGFβ
acturin
RA
TGFβ + acturin
⊖: EB media

(RA: neuronal, Acturin + TGF = muscle)

9. Add 25ml BB media in 50ml Falcon. Add GF's

RA ~~10~~ 25ml (should be 1ml/25ml)

TGFβ 1ml per 1ml media

Acturin 1ml per 1ml media

10. Use 6-well plates "not for tissue culturing".
11. Add 4ml TBS to cells. Resuspend. Centrifuge at 800, 3min
12. Prepare matrigel 1:1 with the right media. 50ml of each.

[Cell Counting: use eastern pipette. Count # cells in 5x5 area. Each square has 16 sq.]

$$\# \text{ cells} \times 10^4 / \text{ml} =$$

counted $200 \times 10^4 = 2 \times 10^6 / \text{ml}$ cells $\times 6 \text{ ml} = 12 \text{ million cells}$

13. Aliquot the cells. Spin down. Add the matrigel media into the cells. Resuspend.