

# EXHIBIT "A"

11/13/01 -11/27/01

EB on Scaffold Expt: Comparison of RA, TGF $\beta$ , and Activin Growth Factors

## METHODS:

### Preparation of Scaffolds:

1. Sterilized 10 polymer scaffolds overnight in 70% ETOH.
2. Soaked scaffolds in 3 changes of PBS, about 5 mins each change. Transferred with sterilized forceps.

### Preparation of EBs:

3. Transfer EBs in media to 50 mL Falcon tube. Pipetted lightly to remove EBs adhered to plate. Washed plate with 5 mL old media.
4. Removed most supernatant. Resuspended in 5 mL new EB media. Transferred to 15 mL Falcon tube.
5. Centrifuged at 800 for 1 minute with brake.
6. Diluted mef trypsin 1:5 (white: 25g porcine tryp/L in 0.9% NaCl) in PBS. Add 2 mL diluted trypsin to the 10mL EB/media mixture. Resuspended by pipetting.
7. Incubated for 5 mins. Added 4mL TNS to cells. Resuspended, and centrifuged at 800 for 3 min. Remove supernatant. Add EB media to aliquot the cells into 2:1 (for scaffolds/EB: EB samples)

### Scaffold Conditions:

| GF              | Stock            | Working   | # scaffolds | EB Plates w/o Scaffold |
|-----------------|------------------|-----------|-------------|------------------------|
| RA              | 7.5mg/mL in DMSO | (1:25000) | 2           | 1                      |
| TGF-B           | 2ug/mL           | (1:1000)  | 2           | 1                      |
| Activin         | 20ug/mL          | (1:1000)  | 2           | 1                      |
| Activin & TGF-B |                  | (1:1000)  | 2           | 1                      |
| Normal media    |                  |           | 2           | 1                      |

8. Prepared 25mL EB media with GFs. (RA 1:25000; TGF-B 1:1000; Act 1:1000)
9. Added 1mL of respective media to each well to be used for scaffolds. Added scaffolds to wells.
10. Prepared matrigel 1:1 with respective media (50 uL matrigel with each media).
11. Aliquoted cells into eppendorfs. Spin down <1000 for 4 min. Remove supernatant carefully. Add 20-25 uL matrigel-media to cells. Mix well.
12. Removed media from scaffolds in wells.
13. Added 20-25 uL cell-matrigel mixture onto each scaffold.
14. Incubate 30 mins to solidify matrigel.
15. Added 4 mL respective media to each well.
16. For remaining EB's aliquot, spin down and add media. Pipette into wells and add respective media.
17. Placed on shaker.

## Ebon Scaffold with GF Expt

Scaffold - PLGA

70/1

1. Sterilize polymer scaffold in ethanol overnight
2. Rinse in PBS. Pour out until only 10mL ethanol. Pour PBS in small wells. Pour all scaffold's ethanol into plate. Transfer into PBS in 6-well plates with scraper.
3. Transfer 3L into 25mL each
4. ~~EB~~ (pyglutathione) Transfer EB media to 50mL Falcon. Pipette lightly to remove EB stuck to plate. Wash with clean media or old media.
5. Remove supernatant until 5mL left. Resuspend in new media (EB media). Transfer to 15mL Falcon. Use same 5mL for all 5 tubes.
6. Repeat again with 5mL S. (centrifuge 800 for 1 min with brake).
7. Dilute Trypsin (Aster 10X): 1 → 5. Add 2mL diluted to the 10mL EB-media mixture. Resuspend by pipetting. ~~Pour into 2 small petri dishes (5mL each).~~
8. Incubate 3min in 15mL Falcon

[Dilute Trypsin 1:5 using PBS → only for EB's]

## Scaffold Conditions - 10

1. RA ~~10~~ 2
2. TGFβ ~~10~~ 2
3. Activin ~~10~~ 2
4. TGFβ + Activin ~~10~~ 2
5. ~~control~~

EB's w/o scaffold

TGFβ  
Activin  
RA  
TGFβ or Activin  
⊕: EB media

(RA: neurons, Activin + TGF = media)

9. Add 25mL EB media in 50mL Falcon. Add GF's  
RA ~~25mL~~ (should be ~~last~~ 25mL)  
TGFβ 1mL per 1mL media  
Activin 1mL per 1mL media
10. Use 6-well plates "not from tissue culture".
11. Add 4mL TNS to cells. Resuspend. Centrifuge at 800, 3min.
12. Prepare material 1:1 with the right media. 50mL of each.

[Cell counting: use eastern pipette, Count # cells in 5x5 area. Each square has 16 sq.]

$$\# \text{ cells} \times 10^4 / \text{mL} =$$

$$\text{counted } 200 \times 10^4 = 2 \times 10^6 / \text{mL cells} \times 6 \text{ mL} = 12 \text{ million cells}$$

13. Aliquot the cells. Spin down. Add the material media into the cells. Resuspend.